



AMERICAN ASSOCIATION OF BIOANALYSTS  
 205 West Levee St. Brownsville, TX 78520  
 (281) 436-5357 - (800) 234-5315 - Fax (713) 781-5008

PARTICIPANT STATISTICS

Immunochemistry

Chemistry Q1 2010

Name	Line No.	Specimen 1			Specimen 2			No. Labs		
		Range & Type	Mean	SD	Range & Type	Mean	SD			
<b>C-Peptide</b>										
<b>Total Population</b>										
Whole Population	1	6.29-11.69	P	8.992	2.38	1.08-2.00	P	1.536	2.03	16



AMERICAN ASSOCIATION OF BIOANALYSTS  
 205 West Levee St. Brownsville, TX 78520  
 (281) 436-5357 - (800) 234-5315 - Fax (713) 781-5008

PARTICIPANT STATISTICS

Immunochemistry

Chemistry Q1 2010

Name	Line No.	Specimen 1		Mean	SD	Specimen 2		Mean	SD	No. Labs
		Range & Type				Range & Type				
<b>Insulin</b>										
<b>Initial Grouping by Reagent and Instrument</b>										
Siemens Immulite 2000, 2500 & Siemens Immulite 2000	1	37.2-69.1	P	53.13	4.7	5.6-10.4	P	8.04	1.0	16
<b>Initial Grouping by Reagent</b>										
Siemens Immulite 2000, 2500	2	37.2-69.1	P	53.13	4.7	5.6-10.4	P	8.04	1.0	16
<b>Initial Grouping by Sensitivity or Principle</b>										
Low recovery methods	3	36.7-68.2	P	52.44	4.4	5.4-10.0	P	7.68	1.1	21
High recovery meth	4	75.7-140.6	P	108.18	9.9	9.4-17.4	P	13.39	1.6	20
<b>Total Population</b>										
Whole Population	5	51.7-96.0	P	73.86	21.4	7.4-13.7	P	10.50	2.4	49



AMERICAN ASSOCIATION OF BIOANALYSTS  
 205 West Levee St. Brownsville, TX 78520  
 (281) 436-5357 - (800) 234-5315 - Fax (713) 781-5008

PARTICIPANT STATISTICS

Immunochemistry

Chemistry Q1 2010

Name	Line No.	Specimen 1			Specimen 2			No. Labs
		Range & Type	Mean	SD	Range & Type	Mean	SD	
<b>Parathyroid Hormone (PTH)</b>								
<b>Initial Grouping by Sensitivity or Principle</b>								
Lower recovery methods	1	124.8-208.0 P	166.40	39.3	21.7-36.2 P	29.00	7.1	29
<b>Total Population</b>								
Whole Population	2	128.7-214.5 P	171.60	25.6	20.7-34.6 P	27.65	4.5	57



AMERICAN ASSOCIATION OF BIOANALYSTS  
 205 West Levee St. Brownsville, TX 78520  
 (281) 436-5357 - (800) 234-5315 - Fax (713) 781-5008

PARTICIPANT STATISTICS

Immunochemistry

Chemistry Q1 2010

Name	Line No.	Specimen 1			Specimen 2			No. Labs
		Range & Type	Mean	SD	Range & Type	Mean	SD	
<b>Vitamin D</b>								
<b>Initial Grouping by Reagent</b>								
DiaSorin LIAISON	1	32.2-77.6 S	54.88	7.6	7.3-21.2 S	14.24	2.3	14
<b>Total Population</b>								
Whole Population	2	0.0-189.5 S	71.19	39.4	0.0-55.2 S	20.03	11.7	26