



**PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS**

205 West Levee St. • Brownsville, TX 78520-5596
800-234-5315 • 281-436-5357 • Fax 713-781-5008

PARTICIPANT STATISTICS

URINE MICROALBUMIN - CREATININE, qualitative

FIRST QUADRIMESTER 2011

Creatinine, semiquantitative

Name	Specimen 1					Specimen 2				
	10 mg/dL	50 mg/dL	100 mg/dL	200 mg/dL	300 mg/dL	10 mg/dL	50 mg/dL	100 mg/dL	200 mg/dL	300 mg/dL
BTNX Rapid Response				1						1
Siemens Clinitek 100				2	6				1	7
Siemens Clinitek 50				1	23				2	22
Siemens Clinitek Status				2	43			1		44
Siemens Micro-Bumintest				1					1	
Total Population	0	0	1	7	75	0	0	2	4	77
Flagging	***	***				***	***			

Microalbumin, semiquantitative or qualitative

Name	Specimen 1								Specimen 2							
	20 mg/L	30 mg/L	50 mg/L	80 mg/L	100 mg/L	150 mg/L	Negative	Positive	0 mg/L	10 mg/L	20 mg/L	30 mg/L	50 mg/L	100 mg/L	Negative	Positive
BTNX Rapid Response				1								1				
Roche Chemstrip Micral			4		13			2		1	7	1	7		2	1
Siemens Clinitek 100						8				8						
Siemens Clinitek 50				5		19				24						
Siemens Clinitek Status				6		43		1		47		2			1	
Total Population	0	0	4	13	15	72	0	4	1	82	8	5	7	0	4	1
Flagging	***	***	***						***				***	***		

Due to a lack of participant consensus, specimen 2 Positive/Negative was not evaluated for this event.

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.