



Whole Blood Glucose LifeScan (WBG or Blood Glucose)

Name	Line No.	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5			No. of Labs					
		Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD						
Initial Grouping by Method																						
LifeScan SureStep	1	140.0 - 210.0	P	174.6	8.4	360.0 - 540.0	P	450.1	30.3	-	-	-	-	-	-	33						
LifeScan SureStep Flexx	2	142.0 - 213.0	P	177.8	14.8	349.0 - 524.0	P	436.3	23.8	67.0 - 101.0	P	83.9	3.8	91.0 - 137.0	P	114.3	5.4	326.0 - 490.0	P	407.9	19.3	361
LifeScan SureStep Pro (P/S)	3	140.0 - 209.0	P	174.4	8.1	342.0 - 513.0	P	427.3	22.3	69.0 - 103.0	P	86	4.4	93.0 - 139.0	P	116	5.5	318.0 - 478.0	P	397.9	20.5	22
Initial Grouping by Sensitivity or Principle																						
LifeScan SureStep (P/S)	4	142.0 - 213.0	P	177.3	14	350.0 - 524.0	P	437	24.7	67.0 - 101.0	P	84.2	4	92.0 - 137.0	P	114.5	5.4	325.0 - 488.0	P	406.3	19.7	416
Total Population																						
Whole Population	5	142.0 - 213.0	P	177.3	14.5	348.0 - 522.0	P	435.2	28.6	67.0 - 101.0	P	84.2	4	92.0 - 137.0	P	114.5	5.4	325.0 - 488.0	P	406.3	19.7	436

Correct responses are defined as those reflecting agreement among 80% or more of all participants.