



Pregnancy

Urine Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alfa Scientific InstantView	10				10		10				10			10	
Beckman ICON 20/25 s/u	32				32		32				32			32	
CardinalHealth hCG Combo	14				14		14				14			13	
Consult Diagnostics	22				22		22				22		1	21	
Quidel QuickVue	18				18		18				18			18	
Sekisui OSOM hCG-Urine	15				15		15				15			15	
Total Population	189	0	0	0	189	0	189	0	0	0	189	0	2	186	0
Flagging		***	***	***		***		***	***	***		***	***		***

Serum Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
CardinalHealth hCG Combo	30				30		30				30			30	
Total Population	89	1	0	0	90	0	88	1	0	0	88	0	0	89	0
Flagging		***	***	***		***		***	***	***		***	***		***

hcg Concentration <2.0 mIU/mL 2815 mIU/mL <2.0 mIU/mL 2815 mIU/mL 2815 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
 Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.