



Urine Drug Screen

Acetaminophen

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Total Population	2	0	2	0
Flagging	***			

Assay Value (mg/dL)

35 mg/dL

Due to a lack of agreement with the target, Specimen 2 was not evaluated this event. The intended result was Present.

Alcohol

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Lin-Zhi	16			16
Syva Emit II, Emit II Plus	12		1	11
Thermo Scientific DRI	18	1	2	17
Total Population	72	1	5	68
Flagging	***		***	

Assay Value (mg/dL)

130 mg/dl

Amphetamines

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	18			18
Bio-Rad TOX/See	11			11
Branan Medical Corporation	17			17
Lin-Zhi	22			22
MedTox Profile V	22			22
Noble One Step	18			17
Syva Emit II, Emit II Plus	15			15
Thermo Scientific DRI	23	1		24
Total Population	296	4	297	2
Flagging	***		***	

Assay Value (ug/mL)

Barbiturates

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	16			16
Bio-Rad TOX/See	11			11
Branan Medical Corporation	18			18
Lin-Zhi	13			13
MedTox Profile V	23			23
Noble One Step	18			17
Syva Emit II, Emit II Plus	12			12
Thermo Scientific DRI	17			17
Total Population	250	0	249	0
Flagging	***		***	

Assay Value (ug/mL)

Benzodiazepines

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere		18		18
Bio-Rad TOX/See		11		11
Branan Medical Corporation	1	17		17
Lin-Zhi		23		22
MedTox Profile V		22		22
Noble One Step		18		17
Syva Emit II, Emit II Plus		11		11
Thermo Scientific DRI	1	20		19
Total Population	4	278	276	5
Flagging	***		***	

Assay Value (ug/mL)

2.9 ug/mL

Buprenorphine

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Branan Medical Corporation	17			17
Lin-Zhi	20			20
MedTox Profile V	13			13
Siemens Viva ser	10			10
Syva Emit II, Emit II Plus	18			18
Total Population	147	0	147	0
Flagging	***		***	

Assay Value (ng/mL)

Cannabinoids

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	18			18
Alfa Scientific InstantView	17			17
Bio-Rad TOX/See	11			11
Branan Medical Corporation	18			18
Lin-Zhi	22			22
MedTox Profile V	23			23
Noble One Step	18			17
Siemens Viva ser	10			10
Syva Emit II, Emit II Plus	32			32
Thermo Scientific DRI	24			24
Total Population	342	1	341	1
Flagging	***		***	

Assay Value (ng/mL)

Cocaine Metabolite - Benzoylcegonine

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	18			18
Alfa Scientific InstantView	17			17
Bio-Rad TOX/See	11			11
Branan Medical Corporation	17			17
CLIAwaived OneScreen	12	1	12	1
Lin-Zhi	25			25
MedTox Profile V	22			22
Noble One Step	18			17
OneStep 3rd party, OEM	10			10
OneStep OEM 3rd party wvd	19			19
Syva Emit II, Emit II Plus	16			16
Thermo Scientific DRI	24			24
Total Population	354	2	353	2
Flagging	***		***	

Assay Value (ug/mL)

Cotinine

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Thermo Scientific DRI		4	4	
Total Population	1	10	11	
Flagging	***			

Assay Value (ng/mL) 650 ng/mL

Fentanyl

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Immunalysis	3		3	
Thermo Scientific DRI	6		6	
Total Population	18	0	18	0
Flagging		***		***

Assay Value (ug/mL) 0.075 ug/mL

K2 Metabolite (UR-144, JWH-073)

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Syva Emit II, Emit II Plus	14		14	
Total Population	40	0	40	0
Flagging		***		***

Assay Value (ng/mL)

Lysergic Acid Diethylamide

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Total Population	3	1	3	1
Flagging				

Assay Value (ug/mL)

*Due to a lack of participant consensus, Specimens 1 and 2 were not evaluated this event. The intended results were negative.

6-Monoacetylmorphine (6-MAM)

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Total Population	5	0	5	0
Flagging		***		***

Assay Value (ug/mL)

MDMA (Ecstasy)

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	11		11	
Corporation	17		17	
Noble One Step	18		17	
Thermo Scientific DRI	12	1	12	1
Total Population	144	1	143	1
Flagging		***		***

Assay Value (ug/mL)

Methamphetamine

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	13		13	
Alfa Scientific InstantView	4		4	
Bio-Rad TOX/See	11		11	
BioView	4		4	
BMC QuickTox	4		4	
Branan Medical Corporation	18		18	
CLIAwaived OneScreen	11	1	11	1
iCup	7		7	
Integrated E-Z Split	3		3	
Lin-Zhi	17		17	
MedTox Profile series	7		7	
MedTox Profile V	19		19	
Noble One Step	18		17	
OneStep 3rd party, OEM	8		8	
OneStep 3rd party, OEM	18		18	
Rapid Response	4		4	
Redwood BioTech One Step	4		4	
Synermed	7		7	
US Diagnostics ProScreen	3		3	
Total Population	212	1	211	1
Flagging		***		***

Assay Value (ug/mL)

Methaqualone

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Thermo Scientific DRI	3			3
Total Population	12	1	2	11
Flagging		***	***	

Assay Value (ug/mL)

2.0 ug/mL

Methadone

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	10		10	
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	17		17	
Lin-Zhi	18		18	
MedTox Profile V	19		19	
Noble One Step	18		17	
Syva Emit II, Emit II Plus	12		12	
Thermo Scientific DRI	15		15	
Total Population	230	0	229	0
Flagging		***		***

Assay Value (ug/mL)

Methadone Metabolite (EDDP)

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Thermo Scientific DRI	4		4	
Total Population	17	0	16	0
Flagging		***		***

Assay Value (ug/mL)

Opiates

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
3rd party cannot determine	5		5	
Alere	17	1	18	
Alere Triage	4		4	
Alfa Scientific InstantView	1	16	16	1
All chromatographic methods	2	1	3	
Beckman Synchron Emit	6		6	
Bio-Rad TOX/See	10	1	11	
BioView	4		4	
BMC QuickTox	3		3	
Branan Medical Corporation	18		18	
Carolina	3		3	
CLIAwaived OneScreen	12	1	12	1
Germaine Aim	7		7	
iCassette	3		3	
iCup	6	1	6	
Integrated E-Z Split	3		3	
Lin-Zhi	25		25	
MedTox Profile series	3	4	7	
MedTox Profile V	2	20	22	
MultiDrug Cup 3rd party, OEM	3		3	
Noble Insta-Screen	4		4	
Noble Multi-CLIN	5		5	
Noble One Step	18		17	
OneStep 3rd party, OEM	10		10	
OneStep OEM 3rd party wvd	19		19	
Rapid Response	3	1	3	1
Redwood BioTech One Step	4		4	
Roche Cobas	9		9	
Siemens Dimension	9		9	
Siemens Viva ser	10		10	
Siemens V-Twin	5		5	
Synermed	9		9	
Syva Emit	3		3	
Syva Emit II, Emit II Plus	31	1	32	
Thermo Scientific DRI	25		25	
US Diagnostics ProScreen	4		4	
Total Population	335	50	380	3
Flagging				***

Assay Value (ug/mL) 200 ng/mL

*Due to a lack of participant consensus, Specimen 1 was not evaluated this event. The intended results was Present.

Propoxyphene

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	8		8	
Beckman Synchron Emit	5		5	
Branan Medical Corp	17		17	
iCassette	3		3	
iCup	3		3	
MedTox Profile series	6		6	
MedTox Profile V	21		21	
Roche Cobas	3		3	
Syva Emit II, Emit II Plus	6		6	
Thermo Scientific DRI	8		8	
Total Population	101	0	101	0
Flagging		***		***

Assay Value (ug/mL)

Tramadol

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Immunalysis	5		5	
Total Population	9	0	9	0
Flagging		***		***

Assay Value (ug/mL)

Zolpidem

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Total Population	5	0	5	0
Flagging		***		***

Assay Value (ug/mL)

Oxycodone

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	14			14
Branan Medical Corporation	17	1	1	17
CLIAwaived OneScreen	8		4	11
Lin-Zhi	22			22
MedTox Profile V	21			21
Noble One Step	18			17
OneStep OEM 3rd party wvd	19			19
Thermo Scientific DRI	21		1	20
Total Population	237	8	5	239
Flagging		***	***	

Assay Value ng/mL

450 ng/mL

Phencyclidine

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	1	11	12	
Alfa Scientific InstantView		4	4	
Beckman Synchron Emit		4	4	
Bio-Rad TOX/See		11	11	
BioView		4	4	
BMC QuickTox	1	2	3	
Branan Medical Corporation	1	16	16	1
Carolina		3	3	
Germaine Aim		4	4	
iCassette		3	3	
iCup		5	5	
Lin-Zhi		8	8	
MedTox Profile series		6	6	
MedTox Profile V		23	23	
MultiDrug Cup 3rd party, OEM	1	2	2	1
Noble Insta-Screen		4	4	
Noble Multi-CLIN		6	6	
Noble One Step		13	13	
Rapid Response		3	3	
Redwood BioTech One Step		4	4	
Roche Cobas		7	7	
Siemens Dimension		9	9	
Synermed		5	5	
Syva Emit II, Emit II Plus		8	8	
Thermo Scientific DRI		15	15	
US Diagnostics ProScreen		3	3	
Total Population	5	229	232	2
Flagging	***			***

Assay Value (ng/mL)

200 ng/mL

Tricyclic Antidepressants

Name	Specimen 1		Specimen 2	
	Absent	Present	Absent	Present
Alere	13		13	
Alere Triage	4		4	
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	17		17	
iCassette	3		3	
iCup	3		3	
MedTox Profile series	5		5	
MedTox Profile V	22		22	
Noble Multi-CLIN	3		3	
Noble One Step	17		17	
Rapid Response	3		3	
Redwood BioTech One Step	3		3	
Thermo Scientific DRI	4		4	
US Diagnostics UScreen	2	1	2	1
Total Population	130	1	130	1
Flagging		***		***

Assay Value (ug/mL)