



**PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS**

205 West Levee St. ♦ Brownsville, TX 78520-5596
800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

SECOND QUADRIMESTER 2011

Serum Pregnancy - Qualitative HCG

Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg
Abbott MEIA / CMIA	3			3			3			3			3		
ABI SureStep Combo	5			5			5			5			5		
Acceva hCG Combo	8			8			8			8			8		
Alfa Sci Instant-View Combo	8			8			8			8			8		
American Diag Pro Combo	6			6			5			6			6		
Beckman ICON 20, 25	15			15			15			15			14		1
Beckman ICON II s/u	6			6			6			5		1	6		
Cardinal Health/SP hCG Comb	18			17		1	18			18			17		1
CDS HCG STAT-PAK ULTRA	1			1			1			1			1		
Consult hCG Combo	3			3			3			3			3		
Fisher Sure-Vue hCG-STAT	26	1		26	1		26	1		26	1		24	1	1
Fisher Sure-Vue s/u hCG	26			26			26			26			26		
Genzyme Contrast	3			3			3			3			3		
Genzyme OSOM hCG-Combo	21			20		1	21			21			20		1
Immunostics Detector Conbi	3			3			3			3			3		
Inverness Clearview	8			8			8			8			8		
Inverness Signify	2			2			2			2			2		
LifeSign Status s/u hCG	2			2			2			2			2		
Mainline Confirms	2			2			2			2			2		
McKesson hCG Combo	7			7			7			7			7		
Meridian ImmunoCard STAT	4			4			4			4			4		
OneStep 3rd party, OEM	7			7			7			7			7		
OneStep Com 3rd part non wv	4			4			4			4			4		
Polymedco Poly stat hCG	4			4			4			4			4		
PSS Select	2			2			2			2			2		
Quidel QuickVue OS Combo	57	1		58			58			57	1		58		
Quidel QuickVue+ OS Combo	25			25			25			25			25		
SA Scientific SAS	7			7			7			7			7		
Select Diag hCG Combo	3			3			3			3			3		
Siemens Immulite	4			4			4			4			4		
Stanbio One-Step	2			2			2			2			2		
Stanbio QuPID PLUS	5			5			5			5			5		
Teco Diag visual	2			2			2			2			2		
Wampole Clearview DUO	5			5			5			5			5		
Total Population	344	3	1	343	2	3	344	2	0	342	4	1	339	3	4
Flagging		*****	*****		*****	*****		*****	*****		*****	*****		*****	*****

hcg Concentration

1728 mIU/mL

469 mIU/mL

1407 mIU/mL

736 mIU/mL

1187 mIU/L

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.

Basic Chemistry