



**PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS**

205 West Levee St. • Brownsville, TX 78520-5596
800-234-5315 • 281-436-5357 • Fax 713-781-5008

PARTICIPANT STATISTICS

SECOND QUADRIMESTER 2011

CARDIAC MARKERS / ISOENZYMES

	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated
CK-2 (CK-MB)										
LifeSign MI	0	1	1	0	1	0	1	0	1	0
Nanogen Cardiac STATus	3	0	3	0	3	0	3	0	3	0
Total Population	3	1	4	0	4	0	4	0	4	0
Flagging		*****		*****		*****		*****		*****

	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated
Myoglobin, qualitative										
LifeSign MI	0	1	1	0	1	0	1	0	1	0
Nanogen Cardiac STATus	0	3	3	0	2	1	3	0	3	0
Total Population	0	5	5	0	4	1	5	0	5	0
Flagging	*****			*****		*****		*****		*****

	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated
Troponin I										
LifeSign MI	0	1	1	0	1		1	0	1	0
Nanogen Cardiac STATus	0	3	3	0	3		3	0	3	0
Total Population	0	5	5	0	4	1	5	0	5	0
Flagging	*****			*****		*****		*****		*****

	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated	Elevated	Unelevated
Troponin T										
Roche Trop T Sensitive	0	3	3	0	3	0	3	0	3	0
Total Population	0	3	3	0	3	0	3	0	3	0
Flagging	*****			*****		*****		*****		*****

Correct responses are defined as those reflecting agreement among 80% or more of all participants. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.