



PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS

205 West Levee St. ♦ Brownsville, TX 78520-5596
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

PARTICIPANT STATISTICS

SECOND QUADRIMESTER 2011

Serum Pregnancy

Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg
Alpha Sci Instant-View Cmbo	6					6	6			6					6
Beckman ICON 20/25 s/u	26					26	26			26					26
Beckman ICON II s/u	4	1				5	5			5					5
CardinalHealth hCG Combo	17	1				18	17	1		17	1				18
Fisher Sure-Vue hCG-STAT	15					15	15			15					15
Fisher Sure-Vue s/u hCG	20					20	20			20					20
Genzyme OSOM hCG-Combo	19					19	19			18		1			19
Henry Schein OneStep	6					6	6			6					6
Inverness Aceava Combo	24					24	24			24					24
McKesson Medi-Lab	9					9	8		1	8		1			9
OneStep 3rd party, OEM	3					3	3			3					3
Polymedco Poly stat hCG	6					6	6			6					6
PSS Consult	3					3	3			3					3
PSS Select Combo	4					4	4			4					4
Quidel Quickvue OS Combo	22					22	22			22					22
Quidel Quickvue+ OS Combo	20					20	20			20					20
SA Scientific SAS	4					4	4			4					4
Siemens Immulite	4					4	4			4					4
Wampole Clearview Combo	3					3	3			3					3
Total Population	490	4	2	0	0	496	492	2	2	490	2	4	0	0	496
Flagging		***	***	***	***			***	***		***	***	***	***	

Urine Pregnancy

Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg
Alpha Sci Instant-View Cmbo	11					11	11			11					11
Beckman ICON 20/25 s/u	45					45	45			44		1	1		44
Beckman ICON II s/u	6					6	6			6					6
Cardinal Health SP	12					12	12			12					12
CardinalHealth hCG Combo	18	2				20	18	2		18	2				20
Diag Sp Pregna-Cert	15	1				16	16			16					16
Fisher Sure-Vue hCG-STAT	34					34	34			33		1	1		33
Fisher Sure-Vue s/u hCG	19					19	18	1		19					19
Fisher Sure-Vue urine hCG	9					9	9			9					9
Genzyme OSOM hCG-Combo	19					19	19			18		1			19
Genzyme OSOM hCG-Urine	16					16	16			16					16
Henry Schein OneStep	11					11	11			11					11
Inverness Aceava Combo	15					15	15			15					15
Inverness Clearview Ultra	6					6	6			6					6
Lifesign Status s/u hCG	4					4	4			4					4
Mainline Confirms hCG	3	1				4	3	1		3	1				4
McKesson Medi-Lab	17					17	17			17					17
METHOD CODE NOT	5					5	5			5					5
OneStep 3rd party, OEM	4					4	4			4					4
PerMaxim OneStep	5					5	5			5					5
Polymedco Poly stat hCG	11					11	11			11					11
PSS Consult	15					15	15			15					15
PSS Select Combo	7					7	7			7					7
Qui QuickVue Urine/One-Step	65				1	64	65			65				1	64
Quidel Quickvue OS Combo	33	1				34	34			34					34
Quidel Quickvue+ OS Combo	35	1				36	36			36					35
Stanbio QuPID	4					4	4			4					4
Stanbio QuPID PLUS	7					7	7			7					7
Stanbio True 20 Plus	12					12	12			12					12
Wampole Clearview Combo	6					6	6			6					6
Total Population	527	8	0	0	2	533	531	4	0	529	3	3	2	3	529
Flagging		***	***	***	***			***	***		***	***	***	***	

hcg Concentration

3030 mIU/mL

<1 mIU/mL

3030 mIU/mL

3030 mIU/mL

<1 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
 Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.