



**PROFICIENCY TESTING SERVICE  
AMERICAN ASSOCIATION OF BIOANALYSTS**

205 West Levee St. ♦ Brownsville, TX 78520-5596  
800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

**SECOND QUADRIMESTER 2016**

**Serum Pregnancy - Qualitative HCG**

Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alere		6			6			6			6			6	
Alere Clearview		4			4			4			4			4	
Alfa Scientific InstantView		11			11			11			11			11	
Beckman ICON 20, 25		8			8			8			8			8	
CardinalHealth hCG Combo		7		1	6		1	6		1	6			7	
Consult hCG Combo		11			11			11			11			11	
Fisher Sure-Vue hCG-STAT		7			7			7			7			7	
Fisher Sure-Vue s/u hCG		14			14			14			14			14	
OneStep 3rd party, OEM		5			5			5			5			5	
Polymedco Poly stat hCG		4			4			4			4			4	
Quidel QuickVue OS Combo		11			11			11			11			11	
Quidel QuickVue+ OS Combo		9			9			9			9			9	
Sekisui OSOM hCG-Combo		32			32			32			32			32	
<b>Total Population</b>	<b>0</b>	<b>181</b>	<b>0</b>	<b>1</b>	<b>180</b>	<b>0</b>	<b>1</b>	<b>180</b>	<b>0</b>	<b>1</b>	<b>180</b>	<b>0</b>	<b>0</b>	<b>181</b>	<b>0</b>
<b>Flagging</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>
<i>hcg Concentration</i>	2095.3 mIU/mL			1220.3 mIU/mL			475.2 mIU/mL			1927.3 mIU/mL			2553.0 mIU/mL		

Correct responses are defined as those reflecting agreement among 80% or more of all participants.  
Unacceptable responses are indicated by "\*\*\*\*" on the Flagging line for each specimen.