



Pregnancy

Urine Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alfa Scientific InstantView		11		11			11			11					11
Beckman ICON 20/25 s/u		16		16			16			16					16
CardinalHealth hCG Combo		12		12			12			12					12
Consult Diagnostics		18		18			18			18			1		17
Fisher Sure-Vue hCG-STAT		13		13			13			13					13
Quidel QuickVue		19		19			18	1		19					19
Sekisui OSOM hCG-Urine		14		14			14			14					14
Total Population	0	189	0	188	0	0	186	1	1	189	0	0	1	188	0
Flagging	***		***		***	***		***	***		***	***	***		***

Serum Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
CardinalHealth hCG Combo		24		24			24			24					24
Total Population	0	87	0	87	0	0	87	0	0	87	0	0	0	87	0
Flagging	***		***		***	***		***	***		***	***	***		***

hcg Concentration 2815 mIU/mL <2.0 mIU/mL <2.0 mIU/mL <2.0 mIU/mL 2815 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
 Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.