



Pregnancy

Urine	Specimen 6			Specimen 7			Specimen 8			Specimen 9			Specimen 10		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alfa Scientific InstantView	11				11			11			11			11	
Consult Diagnostics	19				19			19			19			19	
Quidel QuickVue	19				19			19			19			19	
CardinalHealth hCG Combo	9			1	8		1	8		1	8		1	8	
Sekisui OSOM hCG-Urine	16				16			16			16			16	
Quidel Quickvue+ OS Combo	12				12			12			12			12	
Total Population	169	0	0	1	168	0	3	166	0	1	168	0	6	162	1
Flagging		***	***	***		***	***		***	***		***	***		***

Serum	Specimen 6			Specimen 7			Specimen 8			Specimen 9			Specimen 10		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
CardinalHealth hCG Combo	18				18			18			18			18	
Total Population	84	0	0	0	84	0	0	84	0	0	84	0	1	83	0
Flagging		***	***	***		***	***		***	***		***	***		***

<i>hcg Concentration</i>	<2.0 mIU/mL	55.0 mIU/mL	55.0 mIU/mL	55.0 mIU/mL	55.0 mIU/mL
--------------------------	-------------	-------------	-------------	-------------	-------------

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
 Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.