



Urine Drug Screen

Acetaminophen

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Total Population	4	0	3	1
Flagging		***		

Assay Value (mg/dL)

35 mg/dL

Due to a lack of agreement with the target, Specimen 2 was not evaluated this event. The intended result was Present.

Alcohol

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Lin-Zhi	18		18	
Syva Emit II, Emit II Plus	14		14	
Thermo Scientific DRI	22		22	
Total Population	86	0	86	0
Flagging		***		***

Assay Value (mg/dL)

Amphetamines

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	15	2	1	16
Bio-Rad TOX/See	10			10
Branan Medical Corporation	17			18
Lin-Zhi	22		1	21
MedTox Profile V	23			23
Noble One Step	19			19
Siemens Dimension	10			10
Synermed	12			12
Syva Emit II, Emit II Plus	16			16
Thermo Scientific DRI	26			26
Total Population	303	2	2	304
Flagging		***	***	

Assay Value (ug/mL)

4.5

Barbiturates

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	14	1	14	1
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	19		19	
Lin-Zhi	12		12	
MedTox Profile V	24		24	
Noble One Step	18		18	
Syva Emit II, Emit II Plus	13		13	
Thermo Scientific DRI	19		19	
Total Population	255	1	255	1
Flagging		***		***

Assay Value (ug/mL)

Benzodiazepines

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	16	1	16	1
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	19		19	
Lin-Zhi	21		21	
MedTox Profile V	23		23	
Noble One Step	18		18	
Synermed	13		13	
Syva Emit II, Emit II Plus	13		13	
Thermo Scientific DRI	23		22	
Total Population	287	1	285	1
Flagging		***		***

Assay Value (ug/mL)

Buprenorphine

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Branan Medical Corporation	17	1	18	
Lin-Zhi		19	19	
MedTox Profile V		13	13	
Synermed		12	12	
Thermo Scientific DRI		31	31	
Total Population	25	128	152	1
Flagging	***			***

Assay Value (ng/mL)

155

Cannabinoids

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	16	1	16	1
Alfa Scientific InstantView	18		18	
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	19		19	
Lin-Zhi	24		24	
MedTox Profile V	24		24	
Noble One Step	20		20	
Siemens Dimension	10		10	
Synermed	10		10	
Syva Emit II, Emit II Plus	19		19	
Thermo Scientific DRI	49		49	
Total Population	351	1	351	1
Flagging		***		***

Assay Value (ng/mL)

Cocaine Metabolite - Benzoylcegonine

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	16	1	16	1
Alfa Scientific InstantView	18		17	1
Bio-Rad TOX/See	10		10	
Branan Medical Corporation	18		18	
CLIAwaived OneScreen	10		10	
Lin-Zhi	25		25	
MedTox Profile V	23		23	
Noble One Step	19		19	
OneStep 3rd party, OEM	10		10	
OneStep OEM 3rd party wvd	18		18	
Siemens Dimension	10		10	
Synermed	13		13	
Syva Emit II, Emit II Plus	17		17	
Thermo Scientific DRI	26		25	1
Total Population	364	1	362	3
Flagging		***		***

Assay Value (ug/mL)

Cotinine

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Thermo Scientific DRI	4		4	
Total Population	13	0	13	0
Flagging		***		***

Assay Value (ng/mL)

Fentanyl

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Immunalysis	3			3
Thermo Scientific DRI	9			9
Total Population	22	0	3	19
Flagging		***	***	

Assay Value (ug/mL)

0.075

K2 Metabolite (UR-144, JWH-073)

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Thermo Scientific DRI	21		21	
Total Population	35	1	35	1
Flagging				

Assay Value (ng/mL)

50

25

*Due to a lack of agreement with the intended result, Specimens 6 and 7 were not evaluated.

Lysergic Acid Diethylamide

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Total Population	7	0	5	2
Flagging		***		

Assay Value (pg/mL)

*Due to a lack of participant consensus, Specimen 7 was not evaluated this event. The intended results were negative.

6-Monoacetylmorphine (6-MAM)

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Total Population	2	19	21	0
Flagging	***			

Assay Value (ug/mL)

25

MDMA (Ecstasy)

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	12			12
Corporation	17			18
Noble One Step	18			18
Thermo Scientific DRI	12			12
Total Population	144	1	0	146
Flagging		***	***	

Assay Value (ug/mL)

4

Methamphetamine

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
3rd party cannot determ manuf	7			7
Alere	11	1		12
Alfa Scientific InstantView	5			5
Bio-Rad TOX/See	10			10
BioView	4			4
BMC QuickTox	4		1	3
Branan Medical Corporation	18		1	18
CLIAwaived OneScreen	9	2	1	10
iCup	7			7
Integrated E-Z Split	3			3
iScreen	3		1	2
Lin-Zhi	16			16
MedTox Profile series	8		1	7
MedTox Profile V	20			20
Noble Multi-CLIN	3			3
Noble One Step	18		1	17
OneStep 3rd party, OEM	8		1	7
OneStep OEM 3rd party wvd	17			17
Rapid Response	3			3
Redwood BioTech One Step	5			5
Synermed	11			11
US Diagnostics UScreen	3			3
Total Population	218	3	7	215
Flagging		***	***	

Assay Value (ug/mL)

2

Methaqualone

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Thermo Scientific DRI	4		4	
Total Population	16	0	16	0
Flagging		***		***

Assay Value (ug/mL)

Methadone

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Alere	9	1	9	1
Branan Medical Corporation	18		17	1
Lin-Zhi	17		17	
MedTox Profile V	20		20	
Noble One Step	17		16	1
Synermed	12		12	
Syva Emit II, Emit II Plus	13		13	
Thermo Scientific DRI	16		16	
Total Population	233	1	225	9
Flagging		***		***

Assay Value (ug/mL)

Methadone Metabolite (EDDP)

Name	Specimen 6		Specimen 7	
	Absent	Present	Absent	Present
Thermo Scientific DRI	5		5	
Total Population	24	1	24	1
Flagging		***		***

Assay Value (ug/mL)

Opiates

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
3rd party cannot determ manuf	7		7		
Alere	16	1	17		
Alere Triage	5		5		
Alfa Scientific InstantView	18		18		
Beckman Synchron Emit	6		6		
Bio-Rad TOX/See	10		10		
BioView	4		4		
BMC QuickTox	4		4		
Branan Medical Corporation	18		18		
CLIAwaived OneScreen	11		11		
Germaine Aim	7		7		
iCassette	4		4		
iCup	7		7		
Integrated E-Z Split	4		4		
Lin-Zhi	27		27		
MedTox Profile series	7		7		
MedTox Profile V	23		23		
Noble Insta-Screen	3		3		
Noble Multi-CLIN	5		5		
Noble One Step	19		19		
OneStep 3rd party, OEM	10		10		
OneStep OEM 3rd party wvd	18		18		
Ortho Vitros	3		3		
Rapid Response	3		3		
Redwood BioTech One Step	5		5		
Roche Cobas	9		9		
Siemens Dimension	10		10		
Siemens V-Twin	5		5		
Synermed	13		13		
Syva Emit II, Emit II Plus	19		19		
Thermo Scientific DRI	50		50		
US Diagnostics ProScreen	3		3		
US Diagnostics UScreen	1	2	3		
Total Population	390	3	393	0	
Flagging		***		***	

Assay Value (ug/mL)

Propoxyphene

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
Alere	9		9		
Beckman Synchron Emit	5		5		
Branan Medical Corporation	18		18		
iCassette	3		3		
MedTox Profile series	6		6		
MedTox Profile V	22		22		
Syva Emit II, Emit II Plus	6		6		
Thermo Scientific DRI	9		9		
Total Population	100	0	100	0	
Flagging		***		***	

Assay Value (ug/mL)

Tramadol

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
Immunalysis	4		4		
Total Population	10	0	11	0	
Flagging		***		***	

Assay Value (ug/mL)

Zolpidem

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
Total Population	6	0	6	0	
Flagging		***		***	

Assay Value (ug/mL)

Oxycodone

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
Alere	12		12		
Branan Medical Corporation	18		19		
CLIAwaived OneScreen	9	1	10		
Lin-Zhi	22		22		
MedTox Profile V	22		22		
Noble One Step	18		18		
OneStep OEM 3rd party wvd	18		18		
Synermed	12		12		
Thermo Scientific DRI	19		19		
Total Population	239	1	241		
Flagging		***		***	

Assay Value ng/mL

Phencyclidine

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
3rd party cannot determ manuf	4		4		
Alere	10		10		
Alere Triage	4		4		
Alfa Scientific InstantView	3		3		
Beckman Synchron Emit	4		4		
Bio-Rad TOX/See	10		10		
BioView	4		4		
BMC QuickTox	4		4		
Branan Medical Corporation	18		18		
Germaine Aim	4		4		
iCassette	3		3		
iCup	6		6		
Lin-Zhi	7		7		
MedTox Profile series	7		7		
MedTox Profile V	23		23		
Noble Insta-Screen	3		3		
Noble Multi-CLIN	5		4		
Noble One Step	14		14		
Ortho Vitros	3		3		
Redwood BioTech One Step	5		5		
Roche Cobas	7		7		
Siemens Dimension	9		9		
Synermed	8		8		
Syva Emit II, Emit II Plus	9		9		
Thermo Scientific DRI	16		16		
US Diagnostics UScreen	3		3		
Total Population	234	1	234	0	
Flagging		***		***	

Assay Value (ng/mL)

Tricyclic Antidepressants

Name	Specimen 6		Specimen 7		
	Absent	Present	Absent	Present	
Alere	2	10	12		
Alere Triage		5	5		
Bio-Rad TOX/See		9	9		
Branan Medical Corporation		18	17		
iCassette		3	3		
iCup		4	4		
MedTox Profile series		5	5		
MedTox Profile V		23	23		
Noble Multi-CLIN	2	1	3		
Noble One Step	12	7	19		
Redwood BioTech One Step		3	3		
Thermo Scientific DRI		4	4		
US Diagnostics UScreen		4	4		
Total Population	18	117	134	0	
Flagging	***			***	

Assay Value (ug/mL)

4.5