



Whole Blood Glucose, Lifescan Whole Blood Glucose, Lifescan

Name	Line No.	Specimen 1				Specimen 2				Specimen 3				Specimen 4				Specimen 5				No. of Labs
		Range & Type	Mean	SD		Range & Type	Mean	SD		Range & Type	Mean	SD		Range & Type	Mean	SD		Range & Type	Mean	SD		
Whole Blood Glucose, Lifescan																						
Initial Grouping by Method																						
LifeScan SureStep	1	84.0 - 126.0 P	105.3	6.3	125.0 - 187.0 P	155.7	4.3	-	-	-	-	-	-	-	-	-	-	-	-	-	21	
LifeScan SureStep Flexx	2	85.0 - 128.0 P	106.5	5.8	341.0 - 511.0 P	425.8	20.7	175.0 - 263.0 P	219.3	19.4	62.0 - 93.0 P	77.6	4.5	283.0 - 424.0 P	353.6	28.0					384	
LifeScan SureStep Pro (P/S)	3	88.0 - 132.0 P	109.7	4.4	348.0 - 521.0 P	434.4	21.6	179.0 - 268.0 P	223.4	12.5	65.0 - 97.0 P	80.6	4.5	283.0 - 425.0 P	354.0	19.9					20	
Initial Grouping by Sensitivity or Principle																						
LifeScan SureStep (P/S)	4	85.0 - 128.0 P	106.8	6.0	335.0 - 502.0 P	418.7	49.2	176.0 - 264.0 P	219.8	18.7	62.0 - 94.0 P	77.9	4.6	283.0 - 424.0 P	353.7	27.1					425	
Total Population																						
Whole Population	5	85.0 - 127.0 P	105.9	8.6	338.0 - 506.0 P	422.0	31.3	170.0 - 256.0 P	213.1	23.9	62.0 - 93.0 P	77.7	5.3	270.0 - 406.0 P	338.1	45.5					461	