



**PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS**

205 West Levee St. ♦ Brownsville, TX 78520-5596
800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

THIRD QUADRIMESTER 2011

Serum Pregnancy - Qualitative HCG

Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg
Abbott MEIA / CMIA	3			3			3			2			2		
ABI SureStep Combo	5			5			5			5			5		
Acceva hCG Combo	8			8			8			8			8		
Alfa Sci Instant-View Combo	9			9			9			9			9		
American Diag Pro Combo	6			6			6			6			6		
Beckman Access	2			2			2			2			2		
Beckman ICON 20, 25	16			15		1	16			15		1	16		
Beckman ICON II s/u	6			6			6			6			6		
Cardinal Health/SP hCG Comb	19			19			19			19			19		
Consult hCG Combo	4			4			4			4			4		
Fisher Sure-Vue hCG-STAT	28			28			28			28			28		
Fisher Sure-Vue s/u hCG	30		1	30		1	30		1	30		1	30		1
Genzyme Contrast	2			2			2			2			2		
Genzyme OSOM hCG-Combo	18			17		1	18			18			18		
Immunostics Detector Conbi	3			3			3			3			3		
Inverness Clearview	10			10			10			10			10		
Inverness Signify	2			2			2			2			2		
LifeSign Status s/u hCG	3			3			3			3			3		
McKesson hCG Combo	7		1	7		1	7		1	7		1	7		1
Meridian ImmunoCard STAT	4			4			4			4			4		
NDC Pro Advantage	2			2			2			2			2		
OneStep 3rd party, OEM	5			5			5			5			5		
OneStep Com 3rd part non wv	3			3			3			3			3		
Polymedco Poly stat hCG	5			5			5			5			5		
Quidel QuickVue OS Combo	58			58			58			58			58		
Quidel QuickVue+ OS Combo	27			27			27			27			27		
SA Scientific SAS	6			6			6			6			6		
Select Diag hCG Combo	3			3			3			3			3		
Siemens Advia	2			2			2			2			2		
Siemens Immulite	3			3			3			3			3		
Stanbio One-Step	3			3			3			3			3		
Stanbio QuPID PLUS	5			5			5			5			5		
Teco Diag visual	2			2			2			2			2		
Wampole Clearview DUO	2			2			2			2			2		
Total Population	356	2	1	349	2	7	356	1	1	349	1	7	355	1	1
Flagging		*****	*****		*****	*****		*****	*****		*****	*****		*****	*****

hcg Concentration

211 mIU/mL

1485 mIU/mL

657 mIU/mL

985 mIU/mL

2259 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.

Basic Chemistry