



Pregnancy

Urine Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alfa Scientific InstantView		11		11			11			11			11		
Beckman ICON 20/25 s/u		12		12			12			12			12		
CardinalHealth hCG Combo		12		12			12			12			12		
Consult Diagnostics		21		21			21			21			21		
Fisher Sure-Vue hCG-STAT		15		15			15			15			15		
Fisher Sure-Vue s/u hCG		10		10			10			10			10		
Quidel QuickVue		18		18			18			18			18		
Quidel Quickvue OS Combo		11		11			11			11			11		
Sekisui OSOM hCG-Urine		15		15			15			15			15		
Total Population	0	195	0	194	1	0	194	1	0	0	195	0	1	194	0
Flagging	***		***		***	***		***	***	***		***	***		***

Serum Name	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
CardinalHealth hCG Combo		27		27			27			27			27		
Total Population	0	94	0	93	1	0	93	1	0	0	94	0	0	94	0
Flagging	***		***		***	***		***	***	***		***	***		***

hcg Concentration 2778 mIU/mL <2.0 mIU/mL <2.0 mIU/mL 2778 mIU/mL 2778 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants.
 Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.