



**PROFICIENCY TESTING SERVICE**  
**AMERICAN ASSOCIATION OF BIOANALYSTS**  
 205 West Levee St. Brownsville, TX 78520-5596  
 800-234-5315 • 281-436-5357 • Fax 713-781-5008

THIRD QUADRIMESTER 2019

## Pregnancy - Basic Chem

Name	Specimen 11			Specimen 12			Specimen 13			Specimen 14			Specimen 15		
	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind
Alfa Scientific InstantView		11			11			11			11			11	
Beckman ICON 20, 25		4			4			4			4			4	
CardinalHealth hCG Combo		8		1	7		1	7			8			8	
Consult hCG Combo	1	8			9			9			9			9	
Fisher Sure-Vue s/u hCG		13			13			13			13			13	
OneStep 3rd party, OEM		4			4			4			4			4	
Quidel QuickVue OS Combo		15			15			15			15			15	
Quidel QuickVue+ OS Combo		10			10			10			10			10	
Sekisui OSOM hCG-Combo		7			7			7			7			7	
Siemens Advia		4			4			4			4			4	
<b>Total Population</b>	<b>1</b>	<b>126</b>	<b>0</b>	<b>1</b>	<b>125</b>	<b>0</b>	<b>1</b>	<b>126</b>	<b>0</b>	<b>0</b>	<b>127</b>	<b>0</b>	<b>0</b>	<b>127</b>	<b>0</b>
<b>Flagging</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>

*hcg Concentration*                      1942.8 mIU/mL                      405.1 mIU/mL                      744.9 mIU/mL                      1670.3 mIU/mL                      2581.4 mIU/mL

Correct responses are defined as those reflecting agreement among 80% or more of all participants.  
 Unacceptable responses are indicated by "\*\*\*\*\*" on the Flagging line of each specimen.