



**PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS**
205 West Levee St. • Brownsville, TX 78520-5596
800-234-5315 • 281-436-5357 • Fax 713-781-5008

PARTICIPANT STATISTICS

CHLAMYDIA / GC / STREP GROUP B ANTIGEN

FIRST QUADRIMESTER 2017

Chlamydia	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer
BD ProbeTec ET/Viper		25			25		24	1		25			25		
Cepheid		16			16		16			16			16		
Gen-Probe APTIMA Combo 2		14			14		13	1		14			14		
Quidel QuickVue		10			10		10			10			10		
Total Population	0	69	0	0	69	0	67	2	0	69	0	0	69	0	0
Flagging	***			***				***	***		***	***		***	

Neisseria gonorrhoeae	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer
BD ProbeTec ET/Viper	25			2	23		25			24	1		2	23	
Cepheid	16				16		16			15	1			16	
Gen-Probe APTIMA Combo 2	14			3	10	1	13	1		14				13	
Total Population	59	0	0	5	53	1	58	1	0	57	2	0	2	56	0
Flagging		***	***	***				***	***		***	***	***	***	

Strep Group B	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos
Total Population	8	1	9	0	9	0	9	0	0	9
Flagging		***		***		***		***		***

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.