



PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS
 205 West Levee St. ♦ Brownsville, TX 78520-5596
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

PARTICIPANT STATISTICS

FIRST QUADRIMESTER 2018

C-Reactive Protein, quantitative

Name	Specimen 1					Specimen 2				No. of Labs
	Line No.	Range	Type	Mean	SD	Range	Type	Mean	SD	
Initial Grouping by Method										
Siemens Dimension Flex	1	0.9 - 1.7	S	1.31	0.13	5.5 - 7.6	S	6.58	0.34	30
Ortho Vitros	2	0 - 14.2	S	3.04	3.73	5.3 - 8.0	S	6.64	0.46	15
Beckman Olympus turb	3	0 - 10.3	S	2.02	2.75	0 - 23.4	S	7.62	5.28	20
Abbott Aeroset/Architect	4	0.9 - 1.7	C	1.29	0.08	5.3 - 6.8	S	6.09	0.25	12
Initial Grouping by Sensitivity or Principle										
All direct quantitative meth	5	0 - 7.1	S	1.65	1.81	3.4 - 9.2	s	6.31	0.97	84
All wide range methods	6	0 - 14.7	S	2.78	3.99	0 - 21.5	s	7.28	4.75	28
Total Population										
Whole Population	7	0 - 11.0	S	2.11	2.97	0 - 13.9	s	6.47	2.46	117