



**PROFICIENCY TESTING SERVICE**  
**AMERICAN ASSOCIATION OF BIOANALYSTS**  
 205 West Levee St. ♦ Brownsville, TX 78520-5596  
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

**PARTICIPANT STATISTICS**

**FIRST QUADRIMESTER 2019**

**VIRAL ANTIGEN SCREEN, Waived**

Influenza Type A Method	Specimen 1			Specimen 2		
	Neg	Pos	Pos for A or B	Neg	Pos	Pos for A or B
Alere	6			6		
Alere i	20			20		
BD Veritor	63	1	1	63		2
Binax NOW Flu A&B, wv	3			3		
Consult Diagnostics	8			8		
Henry Schein One Step +	3			3		
LifeSign Status	5			5		
Quidel Quick Vue Infl A+B	7			6		
Quidel Sofia	27			24	2	1
Sekisui OSOM A&B	15			15		
<b>Total Population</b>	<b>163</b>	<b>1</b>	<b>1</b>	<b>159</b>	<b>2</b>	<b>3</b>
<b>Flagging</b>		***	***		***	

Influenza Type B Method	Specimen 1		Specimen 2	
	Neg	Pos	Neg	Pos
Alere	5	1		6
Alere I	20			20
BD Veritor	63	1		64
Binax NOW Flu A&B, wv	3			3
Consult Diagnostics	8			8
Henry Schein One Step +	3			3
LifeSign Status	5			5
Quidel Quick Vue Infl A+B	7			8
Quidel Sofia	26		1	26
Sekisui OSOM A&B	15			15
<b>Total Population</b>	<b>159</b>	<b>2</b>	<b>1</b>	<b>162</b>
<b>Flagging</b>		***	***	

Respiratory Syncytial Virus Method	Specimen 1		Specimen 2	
	Neg	Pos	Neg	Pos
Alere Binax Now	17			17
Alere i	5			5
BD Veritor	9			9
Quidel QuickVue	8			8
Quidel Sofia	5			5
<b>Total Population</b>	<b>46</b>	<b>0</b>	<b>0</b>	<b>46</b>
<b>Flagging</b>		***	***	

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "\*\*\*\*\*" on the Flagging line of each specimen.