



CoaguChek XS, XS Plus, Basic

Name	Line No.	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5			No. of Labs
		Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	
CoaguChek XS Plus Prothrombin Time 2 spec																	
Initial Grouping by Reagent																	
Roche CoaguChek XS*	1	1.0 - 1.4	P 1.2	0.0	2.7 - 3.6	P 3.17	0.19	-	-	-	-	-	-	-	-	14	
Roche CoaguChek XS Plus, Plus/Pro	2	11.9 - 16.1	P 13.98	0.85	31.3 - 42.3	P 36.8	2.15	-	-	-	-	-	-	-	-	16	
Initial Grouping by Sensitivity or Principle																	
All methods	3	8.0 - 10.8	P 9.38	6.17	20.3 - 27.4	P 23.87	16.45	-	-	-	-	-	-	-	-	27	
Total Population																	
Whole Population	4	8.0 - 10.8	P 9.38	6.17	20.3 - 27.4	P 23.87	16.45	-	-	-	-	-	-	-	-	27	
CoaguChek XS, XS Plus INR 2 spec																	
Initial Grouping by Reagent																	
Roche CoaguChek XS	1	0.93 - 1.4	P 1.168	0.054	2.45 - 3.67	P 3.062	0.174	-	-	-	-	-	-	-	-	101	
Roche CoaguChek XS Plus, Plus/Pro	2	0.93 - 1.39	P 1.157	0.063	2.47 - 3.7	P 3.083	0.134	-	-	-	-	-	-	-	-	23	
Initial Grouping by Sensitivity or Principle																	
All methods	3	0.93 - 1.4	P 1.166	0.056	2.45 - 3.68	P 3.066	0.167	-	-	-	-	-	-	-	-	124	
Total Population																	
Whole Population	4	0.93 - 1.4	P 1.166	0.056	2.45 - 3.68	P 3.066	0.167	-	-	-	-	-	-	-	-	124	

*The CoaguChek XS does not report PT in seconds. Laboratories reporting here have reported on the wrong line and will be flagged in future events.