



PARTICIPANT STATISTICS

SECOND QUADRIMESTER 2011

ORAL FLUIDS HIV-1 ANTIBODIES

anti-HIV-1 Screening	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Ind	Pos	Neg	Ind	Pos	Neg	Ind	Pos	Neg	Ind	Pos	Neg	Ind	Pos
Avioq HIV-1			17			17	17					17	17		
Bio-Rad GS HIV 1/2 Plus O			21			21	21					21	21		
In-house validated method			2			2	2					2	2		
OraQuick ADVANCE Rapid			10			10	10					10	10		
OraQuick Rapid HIV-1			5			5									
Total Population	0	0	55	0	0	55	50	0	0	0	0	50	50	0	0
Flagging	****	*****		****	*****			****	*****	****	*****		****	*****	

anti-HIV-1 Confirmation	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Neg, TNP	Pos	Neg	Neg, TNP	Pos	Neg	Neg, TNP	Pos	Neg	Neg, TNP	Pos	Neg	Neg, TNP	Pos
Bio-Rad HIV-1 Western Blot			2			2	1	1				2	1	1	
bioMerieux Oral Vironostika			1			1	1				1	1			
OraSure Western Blot			37			37	20	17				37	20	17	
Total Population	0	3	40	0	3	40	22	21	0	0	3	40	22	21	0
Flagging	****	*****		****	*****				****	*****	*****				****

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.