



PARTICIPANT STATISTICS

SECOND QUADRIMESTER 2016

C. difficile Toxin / Antigen Detection

**C. Difficile Toxin A
(enterotoxin)**

Method	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Negative	Positive	Positive for A or B	Negative	Positive	Positive for A or B	Negative	Positive	Positive for A or B	Negative	Positive	Positive for A or B	Negative	Positive	Positive for A or B
Alere C df Tox A/B I,II	5			5				1	4	5				1	4
Alere Quik Chek Complet	31			31				7	24	30	1		1	6	24
Alere Tox A/B Quik Chek	9			9				2	7	9				2	7
Illumigene	6			6				2	4	6				2	4
Meridian ImmunoCard	10			10				2	8	10				2	8
Remel Xpect	7			7				2	5	7				2	5
Total Population	76	0	0	76	0	0	0	19	57	74	1	1	1	18	57
Flagging		***	***		***	***	***				***	***	***		

**C. Difficile Toxin B
(cytotoxin)**

Method	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
	Negative	Positive	Negative	Positive	Negative	Positive	Negative	Positive	Negative	Positive
Alere Quik Chek Complet	4		4			4	3	1	1	3
Total Population	9	0	9	0	0	8	8	1	1	7
Flagging		***		***	***		***	***		

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.