



**PROFICIENCY TESTING SERVICE**  
**AMERICAN ASSOCIATION OF BIOANALYSTS**  
 205 West Levee St. w Brownsville, TX 78520-5596  
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

**PARTICIPANT STATISTICS**

**VIRAL ANTIGEN SCREEN**

**SECOND QUADRIMESTER 2018**

Influenza Type A	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5		
	Neg	Pos	Pos for A or B	Neg	Pos	Pos for A or B	Neg	Pos	Pos for A or B	Neg	Pos	Pos for A or B	Neg	Pos	Pos for A or B
Alere	3				3		3			3			3		
Alere i	12	3			15		13	2		14	1		13	2	
BD Veritor	28	3	1		14	17	29	2	1	29	1	2	29	2	1
BD Veritor, non-waived	2	1			3		3			2	1		3		
Binax NOW Flu A&B	9				8		9			9			9		
Henry Schein One Step	4				4		4			4			3		
Meridian ImmunoCard STAT	4				3	1	4			4			4		
OraSure	9				9		9			9			9		
Quidel Quick Vue Infl A+B	4				4		4			4			4		
Quidel Sofia Infl A+B	21			3	16		21			21			21		
Remel Xpect Flu A&B	3				3		3			3			3		
Sekisui OSOM A&B	4				4		4			4			4		
<b>Total Population</b>	<b>116</b>	<b>7</b>	<b>1</b>	<b>3</b>	<b>98</b>	<b>19</b>	<b>119</b>	<b>4</b>	<b>1</b>	<b>118</b>	<b>3</b>	<b>3</b>	<b>118</b>	<b>4</b>	<b>1</b>
<b>Flagging</b>		***	***	***				***	***		***			***	***

Influenza Type B	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
Method	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos
Alere i	11	4		15	11	4		15	9	6
BD Veritor	31		1	18	30	1		31	31	
Quidel Sofia Infl A+B	21			19	21			21	19	2
<b>Total Population</b>	<b>117</b>	<b>5</b>	<b>1</b>	<b>105</b>	<b>116</b>	<b>6</b>		<b>122</b>	<b>113</b>	<b>8</b>
<b>Flagging</b>		***	***			***				***

Adenovirus	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
Method	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos
<b>Total Population</b>	<b>0</b>	<b>3</b>	<b>3</b>	<b>0</b>	<b>3</b>	<b>0</b>	<b>3</b>	<b>0</b>	<b>3</b>	<b>0</b>
<b>Flagging</b>	***			***		***		***		***

Respiratory Syncytial Virus	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5	
Method	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos
Alere Binax Now		23	23		23		23		23	
Alere i		6	6		6		6		6	
BD Veritor		3	3		3		3		3	
BD Veritor, non-waived		3	1	2	2	1	3		3	
Binax NOW RSV		6	6		6		6		6	
Meridian Tru RSV		4	4		3	1	4		4	
Quidel QuickVue		9	9		9		9		9	
Quidel Sofia		8	8		8		8		8	
<b>Total Population</b>	<b>0</b>	<b>72</b>	<b>69</b>	<b>3</b>	<b>70</b>	<b>2</b>	<b>72</b>	<b>0</b>	<b>72</b>	<b>0</b>
<b>Flagging</b>	***			***		***		***		***

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "\*\*\*\*\*" on the Flagging line of each specimen.