



**PROFICIENCY TESTING SERVICE  
AMERICAN ASSOCIATION OF BIOANALYSTS**  
205 West Levee St. w Brownsville, TX 78520-5596  
800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

**PARTICIPANT STATISTICS**

**SECOND QUADRIMESTER 2019**

**STREP GROUP A ANTIGEN SCREEN**

Strep A Screen Method	Specimen 6		Specimen 7		Specimen 8		Specimen 9		Specimen 10		No. of Labs
	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	
Alere BioSta Acceava wvd	4			4	4			4	4		4
Alere i	19			19	19			19	19		19
BD Veritor	6			6	6			6	6		6
Cardinal Health SP	4			4	4			4	4		4
Cepheid GeneXpert	7			7	7			7	7		7
Consult Dipstick, waived	7	1		8	8		1	7	8		8
Fisher Sure-Vue, waived	10			10	10			10	10		10
Henry Schein OneStep+, Wv	3			3	3		1	2	3		3
McKesson, waived	8			8	8			8	8		8
Meridian Diag ImmunoCard	3			3	3			3	3		3
Quidel QuickVue +	6			6	6		1	5	6		6
Quidel QuickVue Dipstick	12			12	12			12	12		12
Quidel QuickVue In-Line	32	1	2	30	32	1	1	32	32	1	33
Quidel Sofia	8			8	8			8	8		8
Sekisui OSOM	9			9	9			9	9		9
Sekisui OSOM Ultra	14			14	14			14	14		14
Sekisui OSOM, waived	19			19	19			19	19		19
Sekisui Rapid Strep A	7			7	6	1		7	7		7
<b>Total Population</b>	<b>194</b>	<b>2</b>	<b>2</b>	<b>193</b>	<b>194</b>	<b>2</b>	<b>4</b>	<b>192</b>	<b>195</b>	<b>1</b>	<b>196</b>
<b>Flagging</b>		<b>***</b>	<b>***</b>			<b>***</b>	<b>***</b>			<b>***</b>	

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "\*\*\*\*\*" on the Flagging line of each specimen.