



CoaguChek Prothrombin Time, Basic

Name	Line No.	Specimen 1				Specimen 2				No. of Labs
		Range & Type	Mean	SD	Range & Type	Mean	SD			
Coaguchek Prothrombin Time 2 spec										
Initial Grouping by Method										
Roche CoaguChek S	1	13.2 - 17.9	P	15.54	0.70	20.6 - 27.9	P	24.26	0.97	10
Initial Grouping by Sensitivity or Principle										
All methods	2	13.2 - 17.9	P	15.54	0.70	20.6 - 27.9	P	24.26	0.97	10
Total Population										
Whole Population	3	13.2 - 17.9	P	15.54	0.70	20.6 - 27.9	P	24.26	0.97	10
Coaguchek INR 2 spec										
Initial Grouping by Method										
Roche CoaguChek S	1	1.34 - 2.01	P	1.675	0.360	3.43 - 5.15	P	4.292	0.740	13
Initial Grouping by Sensitivity or Principle										
All methods	2	1.34 - 2.01	P	1.675	0.360	3.34 - 5.01	P	4.173	0.645	13
Total Population										
Whole Population	3	1.34 - 2.01	P	1.675	0.360	3.34 - 5.01	P	4.173	0.645	13