



PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS
 205 West Levee St. • Brownsville, TX 78520-5596
 800-234-5315 • 281-436-5357 • Fax 713-781-5008

THIRD QUADRIMESTER 2016

Rheumatoid Factor, Quantitative

Name	Line No.	Specimen 1				Specimen 2			Specimen 3			Specimen 4			Specimen 5			No. of Labs
		Range & Type	Mean	SD		Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	Range & Type	Mean	SD	
Rheumatoid Factor, quantitative																		
Initial Grouping by Method																		
Roche Tina-Quant	1	5 - 13 C	8.7	3	5 - 13 C	8.7	3	124 - 230 P	177	14.5	5 - 13 C	8.7	3	124 - 230 P	176.7	14.9	6	
Beck Olympus turbidimetric	2	5 - 13 C	9	2	5 - 13 C	9	2	89 - 165 P	126.6	14.7	5 - 13 C	9	2	89 - 165 P	126.6	14.2	5	
Abbott Aeroset/Architect	3	8 - 16 C	12	6	8 - 16 C	12	6	126 - 235 P	180.4	4.8	8 - 16 C	12.2	5.6	126 - 235 P	180.4	6.1	5	
Initial Grouping by Sensitivity or Principle																		
All Roche methods	4	4 - 12 C	7.6	3.8	4 - 12 C	7.5	3.9	118 - 218 P	167.9	22.5	4 - 12 C	7.4	4	120 - 223 P	171.2	24.9	10	
Total Population																		
Whole Population	5	4 - 12 C	8.1	4.4	4 - 12 C	8.1	4.4	128 - 237 P	182.4	65.2	4 - 12 C	8.5	4.4	138 - 256 P	196.8	84.8	32	