



**PROFICIENCY TESTING SERVICE**  
**AMERICAN ASSOCIATION OF BIOANALYSTS**  
 205 West Levee St. ♦ Brownsville, TX 78520-5596  
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

PARTICIPANT STATISTICS

THIRD QUADRIMESTER 2018

**DIAGNOSTIC IMMUNOLOGY**

Infectious Mononucleosis	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5		No. of Labs
	NEG	POS	NEG	POS	NEG	POS	NEG	POS	NEG	POS	
Alere Aceava		33	33		33		33		33		33
Alere Clearview		19	19		19		19		19		19
ASI		14	14		14		14		14		14
Biokit/Fisher Sure-Vue		15	15		15		15		15		15
BTNX		3	3		3		3		3		3
Cardinal Health SP		16	16		16		1		15		16
Consult Diagnostics		5	5		5		5		5		5
Henry Schein One Step +		3	3		3		3		3		3
LifeSign Status Mono		6	6		6		6		6		6
Polymedco Poly Stat		4	4		4		4		4		4
Quidel QuickVue+		16	16		16		16		16		16
Remel Sera-Test/Color Slide		10	10		10		10		10		10
Sekisui OSOM Mono Ser/Plas		14	14		14		14		14		14
Sekisui OSOM Mono Whole BI		11	11		11		11		11		11
Teco Diagnostics		5	5		5		5		5		5
Thermo Fisher Sure-Vue		18	18		18		18		18		18
<b>Total Population</b>	<b>0</b>	<b>208</b>	<b>208</b>	<b>0</b>	<b>208</b>	<b>0</b>	<b>1</b>	<b>207</b>	<b>0</b>	<b>208</b>	<b>208</b>
<b>Flagging</b>	<b>***</b>			<b>***</b>		<b>***</b>	<b>***</b>		<b>***</b>		

Rheumatoid Factor	Specimen 1			Specimen 2			Specimen 3			Specimen 4			Specimen 5			No. of Labs
	NRx	Rx	WRx	NRx	Rx	WRx	NRx	Rx	WRx	NRx	Rx	WRx	NRx	Rx	WRx	
ASI Latex Slide Test	35			33	1	1		35			35			35		35
Biokit	11			11				11			11			11		11
Cardinal Health SP	5			5				5			5			5		5
Fisher Sure-Vue	24	1		25				25			25			25		25
Immunostics immuno/lex-ra	7			7				7			7			7		7
Siemens RapiTex	3			3				3			3			3		3
Stanbio RaPET	10			10				10			10			9	1	10
Teco Diag latex slide test	6			6				6			6			6		6
<b>Total Population</b>	<b>115</b>	<b>1</b>	<b>0</b>	<b>114</b>	<b>1</b>	<b>1</b>	<b>1</b>	<b>115</b>	<b>0</b>	<b>116</b>	<b>0</b>	<b>0</b>	<b>0</b>	<b>115</b>	<b>1</b>	<b>116</b>
<b>Flagging</b>	<b>***</b>		<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>	<b>***</b>

Antistreptolysin O	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5		No. of Labs
	NEG	POS	NEG	POS	NEG	POS	NEG	POS	NEG	POS	
ASI (slide)		24	24			24	24			24	24
Beckman Olympus turb		6	6			6	6			6	6
Fisher Sure-Vue		5	5			5	5			5	5
Immunostics immuno/lex-aso		3	3			3	3			3	3
Roche Tina-Quant		6	6			6	6			6	6
Stanbio (slide)		4	4			4	4			4	4
Teco Diag latex slide test		4	4			4	4			4	4
<b>Total Population</b>	<b>0</b>	<b>55</b>	<b>55</b>	<b>0</b>	<b>0</b>	<b>55</b>	<b>55</b>	<b>0</b>	<b>0</b>	<b>55</b>	<b>55</b>
<b>Flagging</b>	<b>***</b>			<b>***</b>	<b>***</b>			<b>***</b>	<b>***</b>		