



**PROFICIENCY TESTING SERVICE
 AMERICAN ASSOCIATION OF BIOANALYSTS
 205 West Levee St. ♦ Brownsville, TX 78520-5596
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008**

PARTICIPANT STATISTICS

THIRD QUADRIMESTER 2018

HELICOBACTER PYLORI

<i>H. pylori</i> Method	Specimen 1		Specimen 2		No. of Labs
	Neg	Pos	Neg	Pos	
Alere Clearview	16			16	16
Alere ELISA	3			3	3
Alfa Scientific InstantView	19			19	19
BTNX Rapid Response, wv	6			6	6
Cardinal Health SP	7			7	7
Consult Diagnostics	9			9	9
Fisher Sure-Vue IgG	8			8	8
Henry Schein One Step	7			7	7
Micro Detect IgG	4			4	4
Quidel QuickVue gII	22			22	22
Total Population	125	0	0	125	125
Flagging		***	***		

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.