



PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS
 205 West Levee St. w Brownsville, TX 78520-5596
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

PARTICIPANT STATISTICS

THIRD QUADRIMESTER 2018

STREP GROUP A ANTIGEN SCREEN

Strep A Screen Method	Specimen 1		Specimen 2		Specimen 3		Specimen 4		Specimen 5		No. of Labs
	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	
Alere BioSta Acceava wvd	9			9	9			9		9	9
Alere i	7			7	7			7		7	7
BD Veritor	5			5	5			5		5	5
BTNX Rapid Response	3			3	3			3		3	3
Cardinal Health SP	7			7	7			7		7	7
Consult Dipstick, waived	11			11	10	1		11	2	9	11
Fisher Sure-Vue, waived	10			10	10			10		10	10
Henry Schein OneStep+, Wv	4			4	4			4		4	4
McKesson, waived	6			6	6			6		6	6
Meridian Diag ImmunoCard	3			3	3			3		3	3
Quidel QuickVue +	16			16	16			16		16	16
Quidel QuickVue Dipstick	15			15	15			15		15	15
Quidel QuickVue In-Line	42		3	39	41	1	2	40	4	38	42
Quidel Sofia	11			11	10	1		11		11	11
Sekisui OSOM	9			9	9			9		9	9
Sekisui OSOM Ultra	17			17	17			17		17	17
Sekisui OSOM, waived	23			23	23			23		23	23
Sekisui Rapid Strep A	4			4	4			4		4	4
Total Population	218	0	3	215	214	4	2	216	6	212	218
Flagging		***	***			***	***		***		

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.