



PROFICIENCY TESTING SERVICE
AMERICAN ASSOCIATION OF BIOANALYSTS
 205 West Levee St. ♦ Brownsville, TX 78520-5596
 800-234-5315 ♦ 281-436-5357 ♦ Fax 713-781-5008

PARTICIPANT STATISTICS

THIRD QUADRIMESTER 2019

CHLAMYDIA / GC / STREP GROUP B ANTIGEN

Chlamydia	Specimen 11			Specimen 12			Specimen 13			Specimen 14			Specimen 15		
	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer
BD Max	1	2			3		3			2	1		3		
BD ProbeTec ET/Viper		16			16		16			16			16		
Cepheid		24			24		21	1		23	1		19		
Gen-Probe APTIMA Combo		14			14		14			14			14		
Quidel QuickVue		5		1	4		5			5			5		
Total Population	1	63	0	1	63	0	61	1	0	62	2	0	59	0	0
Flagging	***			***				***	***		***			***	***

Neisseria gonorrhoea	Specimen 11			Specimen 12			Specimen 13			Specimen 14			Specimen 15		
	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer	Neg	Pos	Positive for Chlamydia or GC, would refer
BD Max	3				3		3				3			3	
BD ProbeTec ET/Viper	16				16		16				16			16	
Cepheid	24				24		22				24			19	
Gen-Probe APTIMA Combo	14			2	12		14			1	13			14	
Total Population	59	0	0	2	57	0	57	0	0	1	58	0	0	54	0
Flagging		***		***				***	***	***				***	***

Strep Group B	Specimen 11		Specimen 12		Specimen 13		Specimen 14		Specimen 15	
	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos	Neg	Pos
BD Max	3		3		3		3		3	
Cepheid	6		6		6		6		5	
Total Population	11	0	11	0	11	0	0	11	10	0
Flagging		***		***		***	***			***

Correct responses are defined as those reflecting agreement among 80% or more of all participants or referees. Unacceptable responses are indicated by "*****" on the Flagging line of each specimen.